

## CentriPure MINI TRIS Z-50

For rapid purification of proteins or nanoparticles, with simultaneous buffer exchange to TRIS (1 mM, pH 8)



### Instructions for use

#### CentriPure MINI TRIS Z-50

For purification of proteins larger than 25 kD or nanoparticles > 4 nm Ø, with simultaneous buffer exchange to TRIS (1 mM, pH 8).

Rapidly eliminates > 99.999 % of salts, haptens, dyes and other small molecules in less than 5 minutes.

Autoclaved, preservative-free and ready-to-use.

The CentriPure MINI TRIS Z-50 spin column is designed for simultaneous purification and exchange of any buffer to TRIS after a 2 minute centrifugation. Proteins or spheroidal nanoparticles are eluted with minimal dilution into 1 mM TRIS (tris(hydroxymethyl)aminomethane), pH 8.

The protocol provided is optimized for a sample volume of 50 µL. Samples with volumes of 2–100 µL can be processed, but may require further optimization by the end-user.

The gel matrix in the column is Zetadex-50, a size-exclusion gel with an effective pore size of about of 25 kD or about 4 nm Ø. Molecules and particles larger than the pores are excluded from entering the beads, remain in the void volume, pass rapidly through the column, and are eluted free from low molecular weight contaminants.

### Product Parameters


Gel matrix:	Zetadex-50SF beaded dextran
Gel bed volume:	0.5 mL
Storage buffer:	Autoclaved, 1 mM TRIS, pH 8, in sterile water
Sample volume:	2 to 100 µL. Optimized for 50 µL.
Optimal centrifuge conditions:	1000 x g for 2 minutes
Centrifuge type:	Fixed angle or swinging bucket
Removal of fluorescein:	(50 µL 1 mM 5/6-carboxyfluorescein in 0.1 M NaHCO <sub>3</sub> ) > 99.5 %
Removal of fluorescein:	(100 µL 1 mM 5/6-carboxyfluorescein in 0.1 M NaHCO <sub>3</sub> ) > 99 %
Recovery:	Will vary, depending on sample and buffer conditions
Recovery of Dextran Blue:	(50 µL 1 mg/mL Dextran Blue in water) > 70 %
Recovery of Dextran Blue:	(100 µL 1 mg/mL Dextran Blue in water) > 90 %
Size of eluted proteins:	> 25 kD
Size of eluted nanoparticles:	> 4 nm Ø
Storage:	Ambient. DO NOT FREEZE!

#### Materials Provided

- A** CentriPure MINI spin columns (autoclaved and ready-to-use)
- B** Wash Tubes (2 mL)
- C** Sample collection tubes (Blue, 1.5 mL, sterile and RNase/DNase/DNA/pyrogen free)

#### Additional Required Equipment

- Microcentrifuge (capable of generating 1000 x g)
- Variable or fixed pipette (for volumes up to 100 µL)
- Pipette tips
- Microtube rack
- Vortex mixer



If you need assistance or have questions

please call us at +49 30 9489 2201 between 8:00 and 17:00 CET or send an e-mail to [tech@empbiotech.com](mailto:tech@empbiotech.com)


**emp BIOTECH GmbH** · Robert-Rössle-Str. 10 · 13125 Berlin · Germany

**emp BIOTECH GmbH** is an ISO 9001:2015 certified company registration number 011001300789 (TÜV Rheinland)

# CentriPure MINI TRIS Z-50

For rapid purification of proteins larger than 25 kD or nanoparticles >4 nm Ø, with simultaneous buffer exchange to TRIS (1 mM, pH 8)

- robust screw cap
- orientation mark
- ventilation baffles
- porous hydrophilic filter frit
- luer fitting
- snap-off re-usable cap



advanced tapered column design

clicks into wash tube

allows easy removal from centrifuge

## Instructions for use

### 1. Column preparation

- If the columns have been stored cold, allow to warm to room temperature before use (at least 30 minutes).
- Tap gently or briefly vortex to resuspend gel and remove air bubbles.
- Remove CentriPure MINI column from wash tube, loosen the top cap one quarter turn and then snap off the bottom cap.

### 2. Removal of storage buffer

- Click the column back into the wash tube and place assembly in centrifuge. Align the column position upwards using the orientation mark (see Fig. 1).
- Centrifuge at 1000 x g for 2 minutes.\*
- Carefully remove column assembly from centrifuge. Discard wash tube and eluted buffer. Keep the CentriPure MINI column for the next step. Do not disturb gel bed. Use immediately.

### 3. Sample processing

- Carefully apply sample directly to center of gel bed but without touching the gel bed surface (see Fig. 2).
- Place column into a sample collection tube. Insert assembly into the centrifuge. Maintain the upward column alignment (see Fig. 3).
- Centrifuge at 1000 x g for 2 minutes to elute the purified sample.\*
- Remove column assembly from centrifuge. Collect sample and discard column.

\* Calculate the correct rpm speed for 1000 x g by measuring the radius of your centrifuge rotor (center to outer edge of column) using the following formula:

$$\text{rpm} = \sqrt{\frac{1000}{(1.119 \times 10^{-5}) (r)}}$$

Where rpm = revolutions per minute; 1000 represent 1000 x g; and r = radius in cm, measured from center of rotor to outer edge of column.

Example: When the radius, r, is measured to be 7.5 cm, then

$$\text{rpm} = \sqrt{\frac{1000}{(1.119 \times 10^{-5}) (7.5)}} = 3450 \text{ rpm}$$

Figure 1

Orientation mark on column is aligned upwards for storage buffer elution




Figure 2

Proper sample application is made directly onto the center of the gel bed, without touching the gel bed surface





Figure 3

Orientation mark on column is aligned upwards for sample elution



The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final fitness-for-use must be determined by and is the sole responsibility of the end-user.