

Protein and Oligo Labeling

Dyes for protein
and oligonucleotide
labeling



*excellence
made possible*

Protein and Oligo Labeling

Dyes for protein and oligonucleotide labeling



Protein and Oligo Labeling	5
TurboTag™ Labeling Kits	11
Terms and Conditions	14

Protein and Oligo Labeling




emp BIOTECH is a primary manufacturer of a wide variety of standard and specialty dyes and labels.

Fluorophores are essential for use in a vast array of detection methods but have notoriously inconsistent purity profiles. This can negatively affect its extinction coefficient and the sensitivity of the detection method.

emp BIOTECH manufactures high quality NHS esters for post-synthetic modifications for convenient covalent labeling of nucleic acids at the 5'- or 3'-end. Our products are recognized worldwide for their consistent high purity and quality.

emp BIOTECH manufactures dyes and labels under ISO 9001-certified processes, with strict quality controls to ensure extinction coefficients are above specification requirements. This delivers reliable and consistent fluorophore performance, irrespective of the application.

Fluorophore Activated Esters




Amino-functionalized oligonucleotides react with succinimidyl ester activated dyes to form conjugates having stable amide bonds.

TECHNICAL NOTE: Before starting the coupling reaction, the oligonucleotide must be lyophilized once with 0.1 M sodium bicarbonate in order to remove residue ammonia ions!

Product	Order No.	Unit Size
MANT-NHS ester Exc.: 325 nm Em.: 412 nm	AF-0401-D005.0-001	5 mg
5-FAM-NHS ester Exc.: 496 nm Em.: 519 nm	AF-0105-D005.0-001	5 mg
	AF-0105-D020.0-001	20 mg
	AF-0105-D100.0-001	100 mg
6-FAM-NHS ester Exc.: 497 nm Em.: 519 nm	AF-0106-D005.0-001	5 mg
	AF-0106-D020.0-001	20 mg
5-TAMRA-NHS ester Exc.: 550 nm Em.: 575 nm	AF-0107-D005.0-001	5 mg
	AF-0107-D020.0-001	20 mg
6-TAMRA-NHS ester Exc.: 549 nm Em.: 572 nm	AF-0108-D005.0-001	5 mg
	AF-0108-D020.0-001	20 mg
6-ROX-NHS ester Exc.: 587 nm Em.: 599 nm	AF-0110-D005.0-001	5 mg
	AF-0110-D020.0-001	20 mg

Biosensoric








Amino-functionalized oligonucleotides are reacted with succinimidyl ester activated redox labels to form conjugates having stable amide bonds.

TECHNICAL NOTE: Before starting the coupling reaction, the oligonucleotide must be lyophilized once with 0.1 M sodium bicarbonate in order to remove residue ammonia ions!

Product	Order No.	Unit Size
Dicarboxymethylene Blue NHS ester (DCMB-SE)	AF-0201-D005.0-001	5 mg
Monocarboxymethylene Blue NHS ester (MCMB-SE)	AF-0202-D005.0-001	5 mg
Ferrocene carboxylic acid NHS ester (Ferrocene-SE)	AF-0402-D005.0-001	5 mg
Ferrocene-amidopentyl carboxylic acid NHS ester (Ferrocene-C6-SE)	AF-0403-D005.0-001	5 mg
Anthraquinone-2-amidopentyl carboxylic acid NHS ester (AQI-C6-SE)	AF-0405-D005.0-001	5 mg

Additional novel ElectroActive biosensors are currently under development.

Biosensoric

Product Code	Class of Redox Label	Electrode*	Potential Ef (V) vs. Ag/AgCl (1M KCl)	Peak Separation ΔEp (mV) [U = 0.1 V/s]	Stability of Electrochemistry [20 Cycles]	Cyclic Voltamogram Au Electrode (black) Au-dsDNA-Electrode (blue)
AF-0201	Phenothiazine	Au	-0.216	31	stable	
		Au-dsDNA	-0.222	81	stable	
AF-0202	Phenothiazine	Au	-0.198	29	stable	
		Au-dsDNA	-0.201	31	stable	
AF-0402	Metal Complex	Au	+0.292	55	stable	
		Au-dsDNA	+0.292	58	stable	
AF-0403	Metal Complex	Au	+0.390	59	stable	
		Au-dsDNA	+0.390	57	stable	
AF-0405	Quinone	Au	-0.425	68	stable	
		Au-dsDNA	-0.428	79	stable	

* Au electrode – gold wire; Au-dsDNA electrode – gold wire modified with dsDNA (18 base pairs)

For electrochemical applications, **emp BIOTECH** now offers a wide variety of ElectroActive biosensors for labeling of biomolecules, specifically nucleic acids, to form labeled probes equipped with redox-active reporter groups.

Amidite Fluorophores and Labels

Product	Order No.	Unit Size
5-Carboxytetramethylrhodamine DMTr-CE-phosphoramidite (5-TAMRA-DMTr-phosphoramidite)	PF-0102-C100.0-001	100 µmol
	PF-0102-D250.0-001	250 mg
6-Carboxytetramethyl-rhodamine DMTr-CE-phosphoramidite (6-TAMRA-DMTr-phosphoramidite)	PF-0103-C100.0-001	100 µmol
	PF-0103-D250.0-001	250 mg
6-Carboxyfluorescein-dipivaloyl CE-phosphoramidite (6-FAM-phosphoramidite)	PF-0104-C100.0-001	100 µmol
	PF-0104-D250.0-001	250 mg
	PF-0104-E001.0-001	1 g

Labels

Product	Order No.	Unit Size
DAPI aqueous solution, 1 mg/mL	F-0410-M001.0-001	1 mL
Hoechst 33342, 20 mM in water	F-0409-M005.0-001	5 mg
Hoechst 33342, 5 mg/mL in water	F-0409-M020.0-001	20 mg

TurboTag™ Labeling Kits


Precision protein labeling kits

10 min.
protocol



TurboTag™ Labeling Kits by **emp BIOTECH** deliver precision labeling and rapid purification of antibodies, enzymes and other proteins. For controlled labeling and purification of proteins having a molecular weight greater than 25 kDa.

TurboTag labeling kits include the unique emp ZetaSpin Columns for efficient and easy purification



TurboTag™ Labeling Kits

- QUICK:** Pure, labeled protein in just 10 minutes
- EASY:** Simply mix, incubate, and purify
- PRECISE:** Choose and control the degree of labeling of your protein
- CONVENIENT:** Includes ready-to-use ZetaSpin columns

chemistry

Protein with Free Amines + NHS-Ester of the Label → Stable Amide Linkage of the Conjugate

typical workflow

1. buffer protein solution
2. dissolve reactive dye
3. mix and incubate 5 min.
4. load the column
5. centrifuge the column
6. conjugate ready to use

- 3 labeling reactions per kit
- Label up to 15 nmol of protein per reaction using covalent NHS-ester chemistry
- Optimized IgG labeling protocols for 100 µg and 1 mg

Fluorescent Dye Data						
Fluorophore	Excitation (λmax in nm)		Emission (λmax in nm)		as an alternative to	Product Code
MANT	331		426		AlexaFluor™ 350	MK-F0108
DY-405	405		423			MK-D0113
DY-415	418		467			MK-D0114
DY-485XL	485		560			MK-D0143
DY-490	491		515		AlexaFluor™ 488	MK-D0125
DY-495	497		523			MK-D0109
Fluorescein (FAM)	498		522			MK-F0101 MK-F0103*
DY-481XL	515		650			MK-D0103
DY-521XL	523		668			MK-D0126
DY-530	539		561			MK-D0127
DY-555	547		572		AlexaFluor™ 546	MK-D0128
DY-554	551		572		CY3™, AlexaFluor™ 555	MK-D0101
Tetramethylrhodamine (TAMRA)	557		574			MK-T0102
DY-550	558		578			MK-D0111
DY-547P1	559		575			MK-D0112
DY-549P1	560		575			MK-D0116
DY-590	580		599		AlexaFluor™ 568	MK-D0102
X-Rhodamine (ROX)	587		599			MK-R0103
DY-594	594		615		AlexaFluor™ 594	MK-D0129
Texas Red	595		615			MK-R0102
DY-634	635		658			MK-D0107
DY-633	637		657			MK-D0104
Semper Red 647	647		665			MK-F0102
DY-647P1	652		663		CY5™, AlexaFluor™ 647	MK-D0110
DY-648P1	653		672			MK-D0115
DY-649P1	654		672			MK-D0117
DY-652	654		675			MK-D0118
DYQ-661	662		n/a**			MK-D0105
Methylene Blue (DCMB)	667		696			MK-E0150
DY-677	673		694		CY5.5™, AlexaFluor™ 680	MK-D0119
DY-675	680		699			MK-D0106
DY-682	692		709			MK-D0108
DY-701	706		731			MK-D0120
DY-700	707		730			MK-D0130
DY-734	736		759			MK-D0121
DY-750	747		776		CY7™	MK-D0132
DY-752	748		772			MK-D0123
DY-749P1	759		780			MK-D0122
DY-777	770		788			MK-D0124

* FITC
** Quencher. No Stokes shift.

Terms and Conditions



For conducting business with **emp BIOTECH**, please review our general terms and conditions as listed on our website www.empbiotech.com.


emp BIOTECH GmbH
Robert-Rössle-Str. 10
13125 Berlin · Germany
Tel. +49 (0)30 94 89 22 01
info@empbiotech.com

www.empbiotech.com

emp BIOTECH LLC
151 New Jersey 33, Suite 255
Manalapan, NJ 07726 · USA
Tel. +1 (732) 986-9552
info-usa@empbiotech.com

emp BIOTECH is ISO 9001:2015 and 13845:2016 certified.
Registration number 011001300789 (TÜV Rheinland)

More to Discover



Zetadex Gel Filtration
for rapid purification, desalting,
and buffer exchange

Zetadex Gel Filtration Resin
CentriPure Columns
CentriPure Column Arrays
ZetaSpin Columns
ZetaPlate Gel Filtration Plates
ZetaSep FPLC Desalting Columns
ZetaPrep FPLC Desalting Columns



Biomolecule Purification
Solutions for downstream processing

Chromatography Resins for clarified feed streams
– Affinity Chromatography
– Ion Exchange Purification
– Hydrophobic Interaction Chromatography
– Activated Zetarose Solid Phases

Solutions for unclarified feed streams
– SMART Chromatography™

Chromatography Resins for polishing steps
– Size Exclusion Chromatography
– ZetaSep FPLC Desalting Columns
– CentriPure Desalting Columns
– CentriPure Buffer Exchange Columns

