


CentriPure P96

Gel Filtration Column Array

Simultaneously processes 96 samples

Designed specifically for automated systems using gravity or vacuum

Standard ANSI-SBS microplate footprint




Purifies samples
between
150 and 300 μ l


The **CentriPure P96** Column Array is designed for 96 simultaneous purifications in a convenient microplate format.

Sample volumes between 150 and 300 μ l can be purified using either gravity or light vacuum.

Precision packed with **Zetadex-25** ultrapure dextran gel, it is the preferred method for removal of small molecules such as buffer salts, dyes, and haptens from proteins larger than 5 kD.



Elution profile overlay of ovalbumin (1 mg/ml) and free dye (TAMRA, 0,1 μ mol) in a 200 μ l sample volume.



Desalting of protein solution (1 mg albumin (OvA) in 1 ml 0.8 M NaCl), elution with water (200 μ l sample volume)



CentriPure P96

Gel Filtration Column Array

Simultaneously processes 96 samples

Designed specifically for automated systems using gravity or vacuum

Standard ANSI-SBS microplate footprint

1. Column Preparation

- Carefully remove the desired number of cap strips from the top of the array and then remove the entire bottom sealing foil.
- Allow excess column fluid to drain (via gravity) into a suitable waste reservoir. A vacuum of 950 mbar may be used with a manifold to accelerate this process.

2. Column Washing / Equilibration


- Wash each column 4 times (approx. 5 ml total) with either deionized water or buffer (use the same buffer for both equilibration and elution).
- Allow the wash buffer to drain completely between each aliquot. A vacuum of 950 mbar may be used to speed up the washing process.

3. Sample Application

- Load your samples (up to 300 μ l) to each column of the array. Do not use vacuum for sample application. If the sample volume is less than 150 μ l, add enough wash or equilibration buffer so that the combined volume of each sample equals 150 μ l.

4. Elution

- Using the chart below, determine the pre-run and elution volumes specific for your sample size.
- Load the pre-run volume to each column and let it completely enter the gel bed. Do not use vacuum.
- Place a collection plate for sample collection under the array.
- Load the correct elution volume to each column and elute the purified sample by gravity.



Sample volume	Pre-run volume	Elution volume	Protein recovery*	Salt removed
150 μ l	200 μ l	300 μ l	98 %	99.9 %
200 μ l	150 μ l	350 μ l	98 %	99.6 %
250 μ l	100 μ l	400 μ l	98 %	99.6 %
300 μ l	0 μ l	500 μ l	98 %	98.9 %

* determined using 1 mg/ml OvA in 0.8 M NaCl

