

Product Specification Sheet

Product: **CentriPure 25 – Z25M**

Description: Gel filtration column for purification and desalting of proteins larger than 5 kDa or nucleic acids longer than 10 bp/nt.

Processes sample volumes of 2.5 mL.

Product Code: CP-0504

Quality Parameter	Specification
Flow rate <i>2.5 mL water</i>	< 85 seconds
Salt removal <i>Conductivity of 3.5 mL eluate from a 2.5 mL 0.8 M NaCl application</i>	< 1000 µS/cm (1 %)
Dextran Blue recovery yield <i>Yield of Dextran Blue in 3.5 mL eluate from a 2.5 mL solution of 1 mg/mL Dextran Blue application</i>	> 85 %

Other Product Properties

Column bed: Zetadex-25M, swollen in deionized water
 pH Stability: 2 – 13
 Preservative: ProClin® 150
 Storage: Ambient
 Shipping: Ambient

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.