

# Product Specification Sheet



Product: **CentriPure 384 - 100  $\mu$ L well volume**  
**Hydrated gel filtration plates**

Description: For dye terminator removal and/or desalting of sample volumes between 2 and 8  $\mu$ L

Product Code: CP-0102

Quality Parameter	Specification
<b>Matrix Zetadex-50SF</b>	Standard Zetadex-50SF specifications for water regain, sedimentation, salt removal and dextran blue recovery
<b>Uniform Plate Filling</b>	Uniform gel bed height (visual inspection)

## Other Product Properties

Column bed: Zetadex-50SF, swollen in deionized water  
pH Stability: 2 – 13  
Preservative: None  
Storage: 4 to 8°C. DO NOT FREEZE.  
Shipping: Ambient

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.