

Product Specification Sheet

Product: **CentriPure 96 - 1000 μ L well volume**
Hydrated gel filtration plates

Description: For dye terminator removal and/or desalting of sample volumes between 10 and 40 μ L

Product Code: CP-0125

Quality Parameter	Specification
Matrix Zetadex-50SF	Standard Zetadex-50SF specifications for water regain, sedimentation, salt removal and dextran blue recovery
Uniform Plate Filling	Uniform gel bed height (visual inspection)

Other Product Properties

Column bed: Zetadex-50SF, swollen in deionized water
 pH Stability: 2 – 13
 Preservative: None
 Storage: 4 to 8°C. DO NOT FREEZE.
 Shipping: Ambient

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.