

Product Specification Sheet



Product: **CentriPure 96-800SD25**
Hydrated gel filtration plates, 800 µL well volume

Description: For desalting and removal of small molecules.
Sample volumes between 10 and 50 µL

Product Code: CP-0130

Quality Parameter	Specification
Matrix Zetadex-25SF	Standard Zetadex-25SF specifications for water regain, sedimentation, salt removal and dextran blue recovery
Uniform Plate Filling	Uniform gel bed height (visual inspection)

Other Product Properties

Column bed: Zetadex-25SF, swollen in deionized water
pH Stability: 2 – 13
Preservative: None
Storage: 4 to 8°C. DO NOT FREEZE.
Shipping: Ambient

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.