

# Product Specification Sheet



Product: **CentriPure 96 - 1000  $\mu$ L well volume,  
filled with Zetadex-25**  
Hydrated gel filtration plates

Description: For desalting and dye terminator removal

Product Code: CP-0160

Quality Parameter	Specification
<b>Matrix Zetadex-25SF</b>	Standard Zetadex-25SF specifications for water regain, sedimentation, salt removal and dextran blue recovery
<b>Uniform Plate Filling</b>	Uniform gel bed height (visual inspection)

## Other Product Properties

Column bed: Zetadex-25SF, swollen in deionized water  
pH Stability: 2 – 13  
Preservative: None  
Storage: 4 to 8°C. DO NOT FREEZE.  
Shipping: Ambient

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.