

## EZ Dry ULTRA 0.5 g Moisture Trap

For maintaining anhydrous conditions of solvents up to 100 mL.

Contains: 0.5 g activated molecular sieves in an inert highly permeable porous membrane.

<p><b>Product Code:</b> HR-0130</p>	
<p>Moisture Trap with molecular sieves</p>	
<p><b>Quality Parameter</b></p>	<p><b>Specification</b></p>
<p>Mass molecular sieves:  Equilibrium water adsorption test of molecular sieves Solvent Standard: 1600ppm water content, 20°C</p>	<p>0.3 g – 0.7 g  ≤ 20 ppm water content after 24h</p>
<p>Notice:        Storage: Shipping:</p>	<p>The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user. Material Safety Data Sheets are available upon request at <a href="mailto:info@empbiotech.com">info@empbiotech.com</a> This product has been manufactured, tested and assembled under the control of an ISO 9001 registered quality system.</p> <p>Ambient. Ambient.</p> <p>This product has passed all emp Biotech GmbH quality control specifications.</p>