

## EZ Dry ULTRA 2 g Moisture Trap

For maintaining anhydrous conditions of solvents up to 500 mL.

Contains: 2.0 g activated molecular sieves in an inert highly permeable porous membrane.

<p><b>Product Code:</b> HR-0132</p>	
<p>Moisture Trap with molecular sieves</p>	
Quality Parameter	Specification
<p>Mass molecular sieves:</p> <p>Equilibrium water adsorption test of molecular sieves</p> <p>Solvent Standard: 1600ppm water content, 20°C</p>	<p>1.8 g – 2.2 g</p> <p>≤ 20 ppm water content after 24h</p>
<p>Notice:</p> <p>Storage:</p> <p>Shipping:</p>	<p>The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.</p> <p>Material Safety Data Sheets are available upon request at <a href="mailto:info@empbiotech.com">info@empbiotech.com</a></p> <p>This product has been manufactured, tested and assembled under the control of an ISO 9001 registered quality system.</p> <p>Ambient.</p> <p>Ambient.</p> <p>This product has passed all emp Biotech GmbH quality control specifications.</p>