

Product Specification Sheet

DCI Solution

Activator for Chemical RNA and DNA Synthesis

Catalog Number:	NC-0105
Product Name:	DCI Solution
Description:	4,5-Dicyanoimidazole in anhydrous Acetonitrile; 0.25 Molar solution
CAS Number:	Acetonitrile 75-05-8 DCI 1122-28-7

Quality Assurance Data	Specification
Appearance	Clear, no precipitation
DCI Content	0.23 – 0.27 M
Water Content	≤ 30.0 ppm

Application Note

Activator for chemical RNA and DNA synthesis

Shipping: Ambient

Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.