

# Product Specification Sheet

## ETT Solution

0.5 M 5-Ethylthiotetrazole in anhydrous Acetonitrile

Catalog Number: NC-0109

Product Name: ETT Solution

Description: 5-Ethylthiotetrazole in anhydrous Acetonitrile;  
0.50 Molar solution

CAS Number: Acetonitrile 75-05-8  
ETT 89797-68-2

Quality Assurance Data	Specification
Appearance	Clear, no precipitation
Purity by HPLC at 220 nm	≥ 98.0%
Water Content	≤ 20.0 ppm

### Application Note

Activator for chemical RNA and DNA synthesis

**Shipping:** Ambient  
**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**