

Product Specification Sheet

DCI Solution

0.5 M 4,5-Dicyanoimidazole in anhydrous Acetonitrile

Catalog Number: NC-0113

Product Name: DCI Solution

Description: 4,5-Dicyanoimidazole in anhydrous Acetonitrile;
0.5 Molar solution

CAS Number: Acetonitrile 75-05-8
DCI 1122-28-7

Quality Assurance Data	Specification
Appearance	Clear, no precipitation
DCI Content	0.48 – 0.52 M
Water Content	≤ 30.0 ppm

Application Note

Activator for chemical RNA and DNA synthesis

Shipping: Ambient
Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.