

# Product Specification Sheet

## 20% Diethylamine in Acetonitrile (v/v), anhydrous

Catalog Number: NC-0302

Product Name: 20% Diethylamine in Acetonitrile

Description: 20% Diethylamine in anhydrous Acetonitrile, (v/v)

CAS Number: Diethylamine 109-89-7  
Acetonitrile 75-05-8

Quality Assurance Data	Specification
Appearance	Clear Solution
Refractive Index	1,351 – 1,353
Water Content	≤ 50 ppm

### Application Note

Cleavage reagent for nucleic acid synthesis.

**Shipping:** Ambient

**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**