

Product Specification Sheet



Deblock 3% DCA in DCM

Product Code NC-0401

Description 3% Dichloroacetic acid in Dichloromethane (v/v)

CAS Number Dichloromethane: 75-09-2
Dichloroacetic acid: 79-43-6

QUALITY PARAMETER	SPECIFICATION	
Appearance	Clear Solution	
DCA Content (Titration)	2.7 – 3.3	Vol%
Water Content (KF)	≤ 30	ppm

Application Detritylation reagent for nucleic acid synthesis

Storage Ambient

Shipping Ambient

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes.

Final Fitness-For-Use must be determined by, and is the sole responsibility of, the end-user. Material Safety Data Sheets are available upon request at info@empbiotech.com. This product has been manufactured, tested, and assembled under the control of an ISO 9001 registered quality system.