

Product Specification Sheet

Deblocking Reagent

3% Trichloroacetic acid in Dichloromethane (w/v)

Catalog Number: NC-0404

Product Name: Deblocking Reagent

Description: Anhydrous mixture containing 3% Trichloroacetic acid in Dichloromethane, (w/v)

CAS Number: Dichloromethane 75-09-2
Trichloroacetic acid 76-09-3

Quality Assurance Data	Specification
Appearance	Clear Solution
Function Test	positive
Water Content	≤ 120.0 ppm

Application Note

Deblocking/detritylization reagent for nucleic acid synthesis

Shipping: Ambient

Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.