

Product Specification Sheet



Acetonitrile, anhydrous

Catalog Number: NC-0602

Product Name: Acetonitrile, anhydrous

Description: Anhydrous Acetonitrile, water-free

CAS Number: 75-05-8

| Quality Assurance Data | Specification |
|------------------------|----------------|
| Acidity Test | Negative |
| Appearance | Clear Solution |
| Purity | ≥ 99.9% |
| Refractive Index | 1,342 – 1,346 |
| Water Content | ≤ 20.0 ppm |

Application Note

High purity wash solvent, anhydrous

Shipping: Ambient

Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.