

Product Specification Sheet



Acetonitrile

Product Code NC-0602

Description ≤ 20 ppm water
CAS Number 75-05-8

QUALITY PARAMETER	SPECIFICATION
Appearance	Clear colorless liquid
Purity (GC)	≥ 99.9 %
Identity (GC)	Conforms
Water Content (KF)	≤ 20 ppm
Acidity as Acetic Acid (Titration)	≤ 0.05 %
Refractive Index	1.342 – 1.346

Application High purity anhydrous solvent

Storage Ambient

Shipping Ambient

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes.

Final Fitness-For-Use must be determined by, and is the sole responsibility of, the end-user. Material Safety Data Sheets are available upon request at info@empbiotech.com. This product has been manufactured, tested, and assembled under the control of an ISO 9001 registered quality system.