

Product Specification Sheet

Pyridine, anhydrous

Catalog Number: NC-0604

Product Name: Pyridine, anhydrous

Description: Anhydrous Pyridine,
Water-free

CAS Number: 110-86-1

Quality Assurance Data	Specification
Appearance	Clear Solution
Purity	≥ 99.95%
Refractive Index	1,508 – 1,512
Water Content	≤ 30.0 ppm
Content Mass(kg)	n/a
Volume (L)	n/a

Application Note

High purity solvent, anhydrous

Shipping: Ambient

Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.