

# Product Specification Sheet

## Acetonitrile anhydrous, $\leq 10$ ppm

Catalog Number: NC-0609

Product Name: Acetonitrile, anhydrous

Description: Anhydrous Acetonitrile, water-free

CAS Number: 75-05-8

Quality Assurance Data	Specification
Acidity Test	Negative
Appearance	Clear Solution
Purity	$\geq 99.9\%$
Refractive Index	1,342 – 1,346
Water Content	$\leq 10.0$ ppm

### Application Note

High purity wash solvent, anhydrous

**Shipping:** Ambient  
**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**