

# Product Specification Sheet

## Solvent Mix - 60% Pyridine and 40% Acetonitrile (v/v), anhydrous

Catalog Number: NC-0612

Product Name: Solvent Mix

Description: 60% Pyridine and 40% Acetonitrile, anhydrous, (v/v)

CAS Number: Pyridine 110-86-1  
Acetonitrile 75-05-8

Quality Assurance Data	Specification
Appearance	Clear Solution
Content Mass (kg)	n/a
Volume (L)	n/a
Refractive Index	1,436 – 1,448
Water Content	≤ 60.0 ppm

### Application Note

High purity solvent, anhydrous

**Shipping:** Ambient  
**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**