

Product Specification Sheet

Capping A

80% Tetrahydrofuran, 10% 2,6-Lutidine and 10% Acetic anhydride (v/v/v)

Catalog Number: NC-0701

Product Name: Capping A

Description: 80% Tetrahydrofuran,
10% 2,6-Lutidine and
10% Acetic anhydride (v/v/v)

CAS Number:	Tetrahydrofuran	109-99-9
	2,6-Lutidine	108-48-5
	Acetic anhydride	108-24-7

Quality Assurance Data	Specification
Appearance	Clear Solution
Refractive Index	1,412 – 1,416
Water content	≤ 100.0 ppm

Application Note

Capping reagent for nucleic acid synthesis.

Shipping: Ambient

Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.