

Product Specification Sheet

Capping A

90% Tetrahydrofuran and 10% Acetic anhydride (v/v)

Catalog Number: NC-0702

Product Name: Capping A

Description: 90% Tetrahydrofuran and
10% Acetic anhydride (v/v)

CAS Number: Tetrahydrofuran 109-99-9
Acetic anhydride 108-24-7

Quality Assurance Data	Specification
Appearance	Clear
Refractive Index	1,404 – 1,406
Water content	≤ 20.0 ppm

Application Note

Capping reagent for nucleic acid synthesis.

Shipping: Ambient

Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.