

# Product Specification Sheet

## Capping A

80% Acetonitrile and 20% 1-Methylimidazole (v/v)

Catalog Number: NC-0705

Product Name: Capping A

Description: 80% Acetonitrile and  
20% 1-Methylimidazole, (v/v)

CAS Number: Acetonitrile 75-05-8  
1-Methylimidazole 616-47-7

Quality Assurance Data	Specification
Appearance	Clear
Refractive Index	1,373 – 1,377
Water content	≤ 50.0 ppm

### Application Note

Capping reagent for nucleic acid synthesis.

**Shipping:** Ambient

**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**