

Product Specification Sheet

Capping A Ultramid

5% Phenoxyacetic anhydride in 85% Tetrahydrofuran and 10% Pyridine (w/v/v)

Catalog Number: NC-0707

Product Name: Capping A Ultramid

Description: 5% Phenoxyacetic anhydride in
85% Tetrahydrofuran and
10% Pyridine (w/v/v)

CAS Number: Tetrahydrofuran 109-99-9
Pac-anhydride 14316-61-1
Pyridine 110-86-1

| Quality Assurance Data | Specification |
|------------------------|----------------|
| Appearance | Clear Solution |
| Refractive Index | 1,423 – 1,428 |
| Water content | ≤ 50.0 ppm |

Application Note

Capping reagent for nucleic acid synthesis

Shipping: Ambient

Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.