

# Product Specification Sheet

## Capping A

75% Acetonitrile and 25% Acetic anhydride (v/v)

Catalog Number: NC-0708

Product Name: Capping A

Description: 25% Acetic anhydride in 75% Acetonitrile (v/v)

CAS Number: Acetonitrile 75-05-8  
Acetic anhydride 108-24-7

Quality Assurance Data	Specification
Appearance	Clear Solution
Refractive Index	1,355 – 1,359
Water Content	≤ 20.0 ppm

### Application Note

Capping reagent for nucleic acid synthesis

**Shipping:** Ambient  
**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**