

# Product Specification Sheet

## Capping B

84% Tetrahydrofuran and 16% 1-Methylimidazole (v/v)

Catalog Number: NC-0801

Product Name: Capping B

Description: 84% Tetrahydrofuran and  
16% 1-Methylimidazole (v/v)

CAS Number: Tetrahydrofurane 109-99-9  
1-Methylimidazole 616-47-7

Quality Assurance Data	Specification
Appearance	Clear Solution
Refractive Index	1,422 – 1,426
Water content	≤ 100.0 ppm

### Application Note

Capping reagent for nucleic acid synthesis.

**Shipping:** Ambient

**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**