

# Product Specification Sheet

## Capping B

80% Tetrahydrofuran, 10% Pyridine and 10% 1-Methylimidazole (v/v/v)

Catalog Number: NC-0803

Product Name: Capping B

Description: 80% Tetrahydrofuran,  
10% Pyridine and  
10% 1-Methylimidazole (v/v/v)

CAS Number: Tetrahydrofuran 109-99-9  
1-Methylimidazole 616-47-7  
Pyridine 110-86-1

Quality Assurance Data	Specification
Appearance	Clear Solution
Refractive Index	1,427 – 1,430
Water content	≤ 50.0 ppm

### Application Note

Capping reagent for nucleic acid synthesis.

**Shipping:** Ambient

**Storage:** Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**