

# Product Specification Sheet

## Capping B1

40% Acetic anhydride in Acetonitrile (v/v)

Catalog Number: NC-0806

Product Name: Capping B1

Description: 60% Acetonitrile and  
40% Acetic anhydride, (v/v)

CAS Number: Acetonitrile 75-05-8  
Acetic anhydride 108-24-7

Quality Assurance Data	Specification
Appearance	Clear Solution
Refractive Index	1,362 – 1,364
Water Content	≤ 30.0 ppm

### Application Note

Capping reagent for nucleic acid synthesis.

Shipping: Ambient

Storage: Ambient

**For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.**