

Product Specification Sheet

Capping B2

60% 2,6-Lutidine in Acetonitrile (v/v)

Catalog Number: NC-0807

Product Name: Capping B2

Description: 60% 2,6-Lutidine and
40% Acetonitrile (v/v)

CAS Number: 2,6-Lutidine 108-48-5
Acetonitrile 75-05-8

Quality Assurance Data	Specification
Appearance	Clear Solution
Refractive Index	1,435 – 1,439
Water Content	≤ 50.0 ppm

Application Note

Capping reagent for nucleic acid synthesis

Shipping: Ambient
Storage: Ambient

For safety instructions see GHS label! The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.