

Zetarose DEAE FF6

Weak Anion Exchange Agarose Resin.

A cross-linked 6 % beaded agarose IEX resin for protein purification.

Product Code	TM-1205
Quality Parameter	Specification:
Ion exchange capacity	0.08 – 0.16 mmol Cl ⁻ / mL resin
Particle diameter 45 to 160 µm	≥ 80 % in range
Other Product Properties	
Solid phase	Zetarose FlashFlow 6 % highly cross-linked beaded agarose, derivatized with diethylaminoethyl.
Particle diameter	100 µm. Particle range 45 - 160 µm (≥ 80 %).
Application	Anionic exchange resin for purification of biomolecules. Zetarose DEAE FF6 is a cross-linked 6 % agarose bead, derivatized with diethylaminoethyl (DEAE) weak anion ligand.
Velocity properties	Linear velocity ≥ 500 cm/h and operating pressure up to 3 bar.
pH stability	2 – 14 (short term), 3 – 12 (long term). The pH stability of the ion exchange resin will be ultimately dependent on the pH stability of the ligand bound to the resin.
Chemical stability	Zetarose DEAE FF6 is generally tolerant of all commonly used aqueous solutions for protein purification. Avoid oxidizing agents, anionic detergents, and buffers.
Storage buffer	20 % ethanol
Storage	+2 to +30 °C. DO NOT FREEZE!
Miscellaneous	---
Notice	The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.