

ZetaCell DEAE Boost

| | |
|---------------------------------|--|
| Product Code | TM-4208 |
| Quality Parameter | Specification |
| Ion exchange capacity | 0.08 – 0.22 mmol Cl ⁻ / mL resin |
| Particle diameter 140 to 360 µm | ≥ 80 % in range |
| | |
| Notice | The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user. |