

LabRack for all CentriPure Gel Filtration Columns



The **emp LabRack** for **CentriPure** Gel Filtration Columns makes purification easy and convenient.

Its unique design allows use with any column size. Washing, elution and sample collection are all performed smoothly and efficiently. Samples may be collected with either a standard 15 ml Falcon tube, a 1.5 ml microcentrifuge tube or a 5 ml Eppendorf Tube®.


The **emp LabRack** is constructed from sturdy *Dibond* composite material, which provides high stability, resistance to solvents, efficient work flow and increased safety.

easy and
convenient



LabRack

for all CentriPure Gel Filtration Columns




Designed for use with all **CentriPure Gel Filtration** columns, the **LabRack System** consists of two separate racks. The upper rack is for securing the **CentriPure Gel Filtration** column, while the lower rack provides space for sample and waste collection.

Waste and washings are collected using 50 ml Falcon tubes.

Depending on the volume required, samples may be collected using either 15 ml Falcon tubes, 5 ml Eppendorf tubes®, or 1.5 ml microcentrifuge tubes.

easy and convenient



03/2014

